

Datasheet: 303001

Description:	GOAT ANTI RAT IgG
Specificity:	IgG
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoprecipitation	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Ig Fraction - liquid
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.0 mg/ml
Immunogen	Normal rat IgG

External Database Links

UniProt:

P20760	Related reagents
P20759	Related reagents
P20761	Related reagents

[P20762](#) [Related reagents](#)

Entrez Gene:

[679045](#) LOC679045 [Related reagents](#)

[299354](#) Ighg [Related reagents](#)

[362795](#) LOC362795 [Related reagents](#)

RRID AB_618757

Specificity **Goat anti rat IgG antibody** is specific for the gamma chain of rat IgG. It may also react with IgG from other species.
Goat anti rat IgG antibody has been cross-adsorbed against rat IgM to minimize cross-reactivity.

Storage -20°C only (ship +4°C)

Guarantee 18 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10123 available at:
10123: <https://www.bio-rad-antibodies.com/uploads/MSDS/10123.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M352200:190320'

Printed on 09 Feb 2021