

Datasheet: 305001

Description:	GOAT ANTI RAT IgG (H/L)
Specificity:	IgG (H/L)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoprecipitation	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Ig Fraction - liquid
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.0 mg/ml
Immunogen	Rat IgG

External Database Links

UniProt:

P01835	Related reagents
P01836	Related reagents
P20766	Related reagents

[P20767](#) [Related reagents](#)
[P20760](#) [Related reagents](#)
[P20759](#) [Related reagents](#)
[P20761](#) [Related reagents](#)
[P20762](#) [Related reagents](#)

Entrez Gene:

[500180](#) LOC500180 [Related reagents](#)
[299354](#) Ighg [Related reagents](#)
[363828](#) RGD1564318 [Related reagents](#)
[679045](#) LOC679045 [Related reagents](#)
[362795](#) LOC362795 [Related reagents](#)

RRID AB_618759

Specificity **Goat anti Rat IgG antibody** reacts with the heavy and light chains of rat IgG, and the light chains of rat IgA and IgM. Cross reactivity with mouse immunoglobulins has been minimised by cross adsorption with pooled mouse sera and plasmacytoma/hybridoma proteins.

References 1. Liu, F. *et al.* (2015) Rapid and reliable healing of critical size bone defects with genetically modified sheep muscle. [Eur Cell Mater. 30: 118-31.](#)

Storage -20°C only (ship +4°C)

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10123 available at: 10123: <https://www.bio-rad-antibodies.com/uploads/MSDS/10123.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M401993:220718'

Printed on 07 Feb 2023