

Datasheet: 305002

BATCH NUMBER 155565

Description:	GOAT ANTI RAT IgG (H/L):FITC
Specificity:	IgG (H/L)
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.1% Sodium azide		
Approx. Protein Concentrations	1.0 mg/ml		
Immunogen	Rat IgG		
External Database	UniProt:		

Links	P01835	Related reagents
	P01836	Related reagents
	P20766	Related reagents
	P20767	Related reagents
	P20760	Related reagents
	P20759	Related reagents
	P20761	Related reagents
	P20762	Related reagents

Entrez Gene:

500180	LOC500180	Related reagents
299354	Ighg	Related reagents
363828	RGD1564318	Related reagents
679045	LOC679045	Related reagents
362795	LOC362795	Related reagents

RRID AB_616789

Specificity **Goat anti Rat IgG antibody** reacts with the heavy and light chains of rat IgG, and the light chains of rat IgA and IgM. Cross reactivity with mouse immunoglobulins has been minimised by cross adsorption with pooled mouse sera and plasmacytoma/hybridoma proteins.

References 1. Liu, F. *et al.* (2015) Rapid and reliable healing of critical size bone defects with genetically modified sheep muscle. [Eur Cell Mater. 30: 118-31.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/305002>
10040

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M373273:200903'

