

Datasheet: 305009

Description:	GOAT ANTI RAT IgG (H/L):RPE
Specificity:	IgG (H/L)
Format:	RPE
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat		
Product Form	Ig Fraction conjugated to R. Phycoerythrin (RPE) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium azide		
Stabilisers	Stabilising agent		
Approx. Protein Concentrations	0.5 mg/ml		
Immunogen	Rat IgG		
External Database Links	UniProt: P01835 Related reagents		

[P01836](#) [Related reagents](#)
[P20766](#) [Related reagents](#)
[P20767](#) [Related reagents](#)
[P20760](#) [Related reagents](#)
[P20759](#) [Related reagents](#)
[P20761](#) [Related reagents](#)
[P20762](#) [Related reagents](#)

Entrez Gene:

[500180](#) LOC500180 [Related reagents](#)
[299354](#) Ighg [Related reagents](#)
[363828](#) RGD1564318 [Related reagents](#)
[679045](#) LOC679045 [Related reagents](#)
[362795](#) LOC362795 [Related reagents](#)

RRID AB_616871

Specificity **Goat anti Rat IgG antibody** reacts with the heavy and light chains of rat IgG, and the light chains of rat IgA and IgM. Cross reactivity with mouse immunoglobulins has been minimised by cross adsorption with pooled mouse sera and plasmacytoma/hybridoma proteins.

References 1. Liu, F. *et al.* (2015) Rapid and reliable healing of critical size bone defects with genetically modified sheep muscle. [Eur Cell Mater. 30: 118-31.](#)

Storage +4°C DO NOT FREEZE
This product is photosensitive and should be protected from light.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10045 available at:
10045: <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350190:190307'

Printed on 23 Mar 2021