

Datasheet: 402001

**BATCH NUMBER 156858**

<b>Description:</b>	GOAT ANTI RABBIT IgM
<b>Specificity:</b>	IgM
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/2000 - 1/4000
Western Blotting			▪	
Immunoblotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**Target Species** Rabbit

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antisera to rabbit IgM were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** Borate buffered saline

**Preservative Stabilisers** None present

**Approx. Protein Concentrations** IgG concentration 1.0 mg/ml

**Immunogen** Rabbit IgM.

---

**External Database**

**Links**

**UniProt:**

[P03988](#)

[Related reagents](#)

[P04221](#)

[Related reagents](#)

---

**RRID**

AB\_609734

---

**Specificity**

**Goat anti Rabbit IgM antibody** recognizes the heavy chain of rabbit IgM and has been cross-absorbed against rabbit IgG. This product may cross react with IgM from other species.

---

**References**

1. Marques, R.M. *et al.* (2012) Early Inflammatory Response of Young Rabbits Attending Natural Resistance to Calicivirus (RHDV) Infection. [Vet Immunol Immunopathol. 150: 181-8.](#)
  2. Saeed, M.I. *et al.* (2015) Systemic antibody response to nano-size calcium phosphate biocompatible adjuvant adsorbed HEV-71 killed vaccine [Clin Exp Vaccine Res. 4: 88-98.](#)
- 

**Storage**

Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10123 available at:  
<https://www.bio-rad-antibodies.com/SDS/402001>  
10123

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M350196:190307'

Printed on 10 Jul 2024