

Datasheet: 402001

BATCH NUMBER 168357

| | |
|----------------------|----------------------|
| Description: | GOAT ANTI RABBIT IgM |
| Specificity: | IgM |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/2000 - 1/4000 |
| Western Blotting | | | ▪ | |
| Immunoblotting | ▪ | | | |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Rabbit

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to rabbit IgM were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution Borate buffered saline

Preservative Stabilisers None present

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Immunogen Rabbit IgM.

External Database

Links

UniProt:

[P03988](#)

[Related reagents](#)

[P04221](#)

[Related reagents](#)

RRID

AB_609734

Specificity

Goat anti Rabbit IgM antibody recognizes the heavy chain of rabbit IgM and has been cross-absorbed against rabbit IgG. This product may cross react with IgM from other species.

References

1. Marques, R.M. *et al.* (2012) Early Inflammatory Response of Young Rabbits Attending Natural Resistance to Calicivirus (RHDV) Infection. [Vet Immunol Immunopathol. 150: 181-8.](#)
 2. Saeed, M.I. *et al.* (2015) Systemic antibody response to nano-size calcium phosphate biocompatible adjuvant adsorbed HEV-71 killed vaccine [Clin Exp Vaccine Res. 4: 88-98.](#)
-

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/402001>
10123

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South

Tel: +1 800 265 7376

America

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M427953:240301'

Printed on 10 Jul 2024
