

Datasheet: 403005

BATCH NUMBER 161986

| | |
|----------------------|--------------------------|
| Description: | GOAT ANTI RABBIT IgG:HRP |
| Specificity: | IgG |
| Format: | HRP |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/2000 - 1/8000 |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|-----------------------|--|
| Target Species | Rabbit |
| Product Form | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid |

Antiserum Preparation Antisera to rabbit IgG were raised by repeated immunisations of goats with normal rabbit IgG.
The Ig fraction was purified by affinity chromatography on pooled rabbit IgG covalently linked to agarose.
The preparation has been cross absorbed against rabbit IgM to remove any potential cross reactivity.

| | |
|------------------------|---------------------------|
| Buffer Solution | Phosphate buffered saline |
|------------------------|---------------------------|

| | |
|---------------------------------|--------------|
| Preservative Stabilisers | 50% glycerol |
|---------------------------------|--------------|

Immunogen Purified rabbit IgG from normal rabbit serum

External Database

Links

UniProt:

[P01870](#)

[Related reagents](#)

Entrez Gene:

[100009097](#)

LOC100009097

[Related reagents](#)

RRID

AB_609686

Specificity

Goat anti Rabbit IgG antibody recognizes the heavy chain of rabbit IgG.

Cross-reactivity against rabbit IgM is reduced by solid phase adsorption. This antibody may react with IgG from other species.

References

1. Wegmann, F. *et al.* (2011) A Novel Strategy for Inducing Enhanced Mucosal HIV-1 Antibody Responses in an Anti-Inflammatory Environment [PLoS One. 6\(1\):e15861.](#)
2. Radomirovic, M. *et al.* (2022) Phycocyanobilin-modified β -lactoglobulin exhibits increased antioxidant properties and stability to digestion and heating [Food Hydrocolloids. 123: 107169.](#)

Storage

Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10106 available at:
<https://www.bio-rad-antibodies.com/SDS/403005>
10106

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350201:190307'

