

Datasheet: 4050-0184

Description:	NATIVE HUMAN ELASTASE (NEUTROPHIL)
Name:	ELASTASE
Other names:	ELA2
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human neutrophils - liquid
Buffer Solution	50mM sodium acetate, pH 5.5 150mM sodium chloride
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.0 mg/ml
External Database Links	UniProt: P08246 Related reagents

Entrez Gene:

[1991](#) ELANE [Related reagents](#)

Synonyms	ELA2
-----------------	------

Product Information	Human elastase is a 28.5 kDa serine protease that has been purified from human neutrophils. Elastase is also present in the pancreas, leucocytes and serum.
----------------------------	---

Protein Molecular Weight	28.5 kDa
---------------------------------	----------

Activity	20-22 U/mg protein 1 unit hydrolyzes one μ mole per minute of MeO-Suc-Ala-Ala-Pro-Val-pNA at 25°C, in 100mM Tris-Hcl, pH 7.5 with 500mM NaCl.
-----------------	--

Purity	>95% by SDS PAGE
---------------	------------------

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10168 available at: https://www.bio-rad-antibodies.com/SDS/4050-0184 10168 Donor material tested and found negative for HBsAg and antibodies to HIV1, HIV2, HCV and HBcAg. As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M362530:200514'

Printed on 25 Mar 2023