

Datasheet: 405002

BATCH NUMBER 156264

Description:	GOAT ANTI RABBIT IgG (H/L):FITC
Specificity:	IgG (H/L)
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/10 - 1/100
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Rabbit

Species Cross Reactivity

Does not react with: Mouse, Human

Product Form

Purified Ig conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid

Max Ex/Em

Fluorophore	Excitation Max (nm)	Emission Max (nm)
FITC	490	525

Antiserum Preparation Antisera to rabbit IgG (H+L) were raised by repeated immunisations of goats with normal rabbit IgG. Purified IgG was prepared from whole serum by affinity chromatography on pooled rabbit IgG covalently linked to agarose. The preparation has been cross absorbed against human and mouse sera and purified paraproteins to minimise potential cross

reactivity.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.1% Sodium Azide (NaN₃)

Approx. Protein Concentrations Ig concentration 1.0 mg/ml

Immunogen Normal rabbit IgG.

External Database Links

UniProt:

[P01839](#) [Related reagents](#)

[P01841](#) [Related reagents](#)

[P01838](#) [Related reagents](#)

[P01847](#) [Related reagents](#)

[P01870](#) [Related reagents](#)

[P01840](#) [Related reagents](#)

Entrez Gene:

[100422811](#) IGKC2 [Related reagents](#)

[100009097](#) LOC100009097 [Related reagents](#)

RRID AB_609712

Specificity **Goat anti Rabbit IgG (H/L) antibody** recognises rabbit IgG heavy and light chains and the light chains of rabbit IgM and IgA as demonstrated by ELISA and flow cytometry.

Flow Cytometry Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.

References 1. Louderback, K.M. *et al.* (2006) Subbarrel patterns of thalamocortical innervation in rat somatosensory cortical barrels: Organization and postnatal development. [J Comp Neurol. 497: 32-41.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/405002>
10040

North & South Tel: +1 800 265 7376

Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M358101:190911'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)