

Datasheet: 4440-0456

Description:	NATIVE HUMAN FIBRIN DEGRADATION PRODUCT E
Name:	FIBRIN DEGRADATION PRODUCT E
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human serum - lyophilized
Reconstitution	Reconstitute with 1.0 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified protein prepared through salt fractionation, gel filtration and ion exchange chromatography.
Buffer Solution	Glycine buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 0.1 mg/ml after reconstitution
External Database Links	UniProt:

[P02679](#) [Related reagents](#)

[P02675](#) [Related reagents](#)

[P02671](#) [Related reagents](#)

Entrez Gene:

[2266](#) FGG [Related reagents](#)

[2244](#) FGB [Related reagents](#)

[2243](#) FGA [Related reagents](#)

Protein Molecular Weight

50 kDa

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10346 available at:
<https://www.bio-rad-antibodies.com/SDS/4440-0456>
10346

Donor material tested and found negative or non-reactive for anti-HIV1/2, HIV1 antigen, HBsAg, syphilis, anti-HCV, anti-HBcAg and anti-HTLV1 & 2.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.

Regulatory

For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M426262:231207'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)