

Datasheet: 4440-0456

Description:	NATIVE HUMAN FIBRIN DEGRADATION PRODUCT E
Name:	FIBRIN DEGRADATION PRODUCT E
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human serum - lyophilized
Reconstitution	Reconstitute with 1.0 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified protein prepared through salt fractionation, gel filtration and ion exchange chromatography.
Buffer Solution	Glycine buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 0.1 mg/ml after reconstitution

External Database Links

UniProt:

P02679	Related reagents
P02675	Related reagents
P02671	Related reagents

Entrez Gene:

2266	FGG	Related reagents
2244	FGB	Related reagents

Protein Molecular Weight

50 kDa

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended.

Shelf Life

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10346 available at:
10346: <https://www.bio-rad-antibodies.com/uploads/MSDS/10346.pdf>

Donor material tested and found negative or non-reactive for anti-HIV1/2, HIV1 antigen, HBsAg, syphilis, anti-HCV, anti-HBcAg and anti-HTLV1 & 2.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.

Regulatory

For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M333204:181012'

Printed on 06 Dec 2018