

Datasheet: 4440-8604

BATCH NUMBER 150729

Description:	NATIVE HUMAN FIBRINOGEN
Name:	FIBRINOGEN
Format:	Purified
Product Type:	Purified Protein
Quantity:	25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human plasma - liquid
Preparation	Protein purified by precipitation
Buffer Solution	Citrate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.25 mg/ml fibrinogen

External Database Links

UniProt:

P02671	Related reagents
P02675	Related reagents
P02679	Related reagents

Entrez Gene:[2243](#) FGA [Related reagents](#)[2244](#) FGB [Related reagents](#)[2266](#) FGG [Related reagents](#)

Product Information **Native Human fibrinogen** is a 340kDa glycoprotein involved in cell adhesion and migration processes including embryogenesis, metastasis, wound healing, blood coagulation and host immunity. It consists of a dimer, each part made up of three disulphide linked polypeptides known as the alpha, beta and gamma chains.

Protein Molecular Weight 340 kDa

Purity 35-65% protein, of which $\geq 90\%$ clottable protein.
Virtually free from plasminogen and plasmin.

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10175 available at:
<https://www.bio-rad-antibodies.com/SDS/4440-8604>
10175

Donor material tested and found negative for HIV1, HIV2, HCV, ALT, Syphilis, HBsAg and antibodies to HCV.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373490:200915'

Printed on 18 Jan 2024