

Datasheet: 4600-1704

BATCH NUMBER 150503

Description:	NATIVE ESCHERICHIA COLI BETA GALACTOSIDASE
Name:	BETA GALACTOSIDASE
Format:	Lyophilized
Product Type:	Purified Protein
Quantity:	5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA			▪	
Western Blotting			▪	
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified protein from <i>Escherichia coli</i> - lyophilized
Reconstitution	Reconstitute with 1.0 ml sterile distilled water or aqueous buffer of choice. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	Contains Sucrose
Approx. Protein Concentrations	Protein content is approximately 20-25% of total weight
External Database Links	UniProt:

Product Information *Escherichia coli* Beta-Galactosidase is an inducible tetrameric enzyme coded by the lac Z gene of the lac operon that is often used as a reporter to assess the efficiency of transfection. It is a metalloenzyme that splits lactose into glucose and galactose. It hydrolyzes terminal, non-reducing beta-D-galactose residues in beta-D-galactosides.

Protein Molecular Weight 540 kDa

Activity >600 U/mg protein
One unit hydrolyses 1 µmole of p-nitrophenol-b-D-galactoside to p-nitrophenol per minute at 37°C.
For batch-specific activity, please inquire.

Storage Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10152 available at:
<https://www.bio-rad-antibodies.com/SDS/4600-1704>
10152

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363127:200528'

Printed on 18 Jan 2024