

Datasheet: 4750-0254

**BATCH NUMBER 310812**

<b>Description:</b>	SHEEP ANTI HUMAN GROWTH HORMONE
<b>Specificity:</b>	GROWTH HORMONE
<b>Other names:</b>	SOMATOTROPIN
<b>Format:</b>	Serum
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.25 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Radioimmunoassays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Serum - liquid
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Immunogen</b>	Purified human growth hormone.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01241</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">2688</a>    GH1    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_618597

<b>Specificity</b>	<b>Sheep anti Human growth hormone antibody</b> recognizes human growth hormone. Growth hormone is a small peptide hormone released by the anterior pituitary under the control of the hypothalamus. It acts to control body growth and metabolic activity.
<b>References</b>	1. Rowe, M. <i>et al.</i> (1994) Upregulation of bcl-2 by the Epstein-Barr virus latent membrane protein LMP1: a B-cell-specific response that is delayed relative to NF-kappa B activation and to induction of cell surface markers. <a href="#">J Virol. 68: 5602-12.</a>
<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10081 available at: <a href="https://www.bio-rad-antibodies.com/SDS/4750-0254">https://www.bio-rad-antibodies.com/SDS/4750-0254</a> 10081
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M363152:200528'

**Printed on 15 Mar 2024**