

Datasheet: 5030-9954

BATCH NUMBER 158567

Description:	NATIVE BOVINE HYALURONIDASE
Name:	HYALURONIDASE
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.5 g

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
-----------------------	--------

Product Form	Purified enzyme from bovine testes - liquid
---------------------	---

Preservative Stabilisers	None present
---------------------------------	--------------

Approx. Protein Concentrations	Total protein concentration 0.2 g/ml
---------------------------------------	--------------------------------------

External Database Links

UniProt:

[Q8SQG8](#) [Related reagents](#)

[Q5E985](#) [Related reagents](#)

Entrez Gene:

[281838](#) HYAL2 [Related reagents](#)

[515397](#) HYAL1 [Related reagents](#)

Product Information	Hyaluronidase catalyzes the depolymerization of mucopolysaccharides, hyaluronic acid,
----------------------------	---

and the chondroitin sulfates A and C. The enzyme is widely distributed in animal tissues but is found in great concentrations in the bovine and ovine testes.

Protein Molecular Weight

55 kDa

Activity

315 U/mg protein

1 unit is defined as the amount of enzyme that liberates one micromole of N-acetylglucosamine per minute at 37°C and pH 4.0.

EC 3.2.1.36

Storage

Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10078 available at:

10078: <https://www.bio-rad-antibodies.com/uploads/MSDS/10078.pdf>

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M362533:200514'

Printed on 07 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)