

Datasheet: 5204-2504

**BATCH NUMBER 154073**

<b>Description:</b>	GOAT ANTI RAT IgG (H/L):HRP
<b>Specificity:</b>	IgG (H/L)
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/500 - 1/1500
ELISA	▪			1/50000 - 1/100000
Western Blotting	▪			1/2000 - 1/10000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Rat
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilised
<b>Reconstitution</b>	Reconstitute with 1.0 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Gentamicin Sulphate 1% Bovine Serum Albumin
<b>Immunogen</b>	Rat IgG, whole molecule.

**External Database  
Links**

**UniProt:**

<a href="#">P01836</a>	<a href="#">Related reagents</a>
<a href="#">P20766</a>	<a href="#">Related reagents</a>
<a href="#">P20760</a>	<a href="#">Related reagents</a>
<a href="#">P20759</a>	<a href="#">Related reagents</a>
<a href="#">P20761</a>	<a href="#">Related reagents</a>
<a href="#">P01835</a>	<a href="#">Related reagents</a>
<a href="#">P20767</a>	<a href="#">Related reagents</a>
<a href="#">P20762</a>	<a href="#">Related reagents</a>

**Entrez Gene:**

<a href="#">363828</a>	RGD1564318	<a href="#">Related reagents</a>
<a href="#">299354</a>	Ighg	<a href="#">Related reagents</a>
<a href="#">500180</a>	LOC500180	<a href="#">Related reagents</a>
<a href="#">679045</a>	LOC679045	<a href="#">Related reagents</a>
<a href="#">362795</a>	LOC362795	<a href="#">Related reagents</a>

---

**RRID** AB\_619913

---

**Specificity** **Horseradish peroxidase conjugated Goat anti Rat IgG antibody** recognizes Rat IgG heavy and light chains. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-peroxidase, anti-goat serum, rat IgG and rat serum.

---

**References** 1. Sanchez de Groot, N. *et al.* (2015) Proteome response at the edge of protein aggregation. [Open Biol. 5. pii: 140221.](#)

---

**Storage** Prior to reconstitution store at +4°C.  
After reconstitution store at +4°C or at -20°C if preferred.  
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10482 available at: <https://www.bio-rad-antibodies.com/SDS/5204-2504>  
10482

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)  
[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)  
[TMB SIGNAL+ \(BUF054A\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M375153:201217'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)