

Datasheet: 5213-2504

**BATCH NUMBER 171298**

<b>Description:</b>	GOAT ANTI RABBIT IgG(Fc):HRP
<b>Specificity:</b>	IgG (Fc)
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/1,000 - 1/5,000
ELISA	▪			1/110,000
Western Blotting	▪			1/2,000 - 1/20,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Ig Fraction conjugated to Horseradish Peroxidase (HRP) - lyophilized
<b>Reconstitution</b>	Use 1.0 ml of sterile DI water. For extended storage reconstitute with 50% glycerol instead of water and then aliquot and store at -20°C. Avoid freeze thaw cycles.
<b>Preparation</b>	Purified Ig prepared by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Gentamicin sulphate 10 mg/ml Bovine Serum Albumin
<b>Immunogen</b>	Native rabbit IgG.

**External Database  
Links**

**UniProt:**

[P01870](#)   [Related reagents](#)

**Entrez Gene:**

[100009097](#)   LOC100009097   [Related reagents](#)

---

**RRID**   AB\_619907

---

**Specificity**   **Goat anti Rabbit IgG (Fc) antibody** recognizes the Fc portion of rabbit IgG. Assay by IEP resulted in a single precipitin arc against anti-peroxidase, anti-goat serum, rabbit IgG, rabbit IgG Fc and rabbit serum. No reaction was observed against rabbit IgG F(ab')<sub>2</sub>.

---

**Further Reading**   1. Farr, A.G. & Nakane, P.K. (1981) Immunohistochemistry with enzyme labeled antibodies: a brief review. [J Immunol Methods. 47 \(2\): 129-44.](#)

---

**Storage**   This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. After reconstitution it is recommended to aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**   12 months from date of despatch

---

**Health And Safety Information**   Material Safety Datasheet documentation #10482 available at: <https://www.bio-rad-antibodies.com/SDS/5213-2504>

---

**Regulatory**   For research purposes only

---

## Related Products

### Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M427987:240301'

Printed on 29 Jan 2026