

## Datasheet: 5268-6010

|                      |                                    |
|----------------------|------------------------------------|
| <b>Description:</b>  | MOUSE ANTI HUMAN KAPPA LIGHT CHAIN |
| <b>Specificity:</b>  | KAPPA LIGHT CHAIN                  |
| <b>Format:</b>       | Purified                           |
| <b>Product Type:</b> | Monoclonal Antibody                |
| <b>Clone:</b>        | 4G7                                |
| <b>Isotype:</b>      | IgG2a                              |
| <b>Quantity:</b>     | 0.5 mg                             |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                            | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Immunohistology - Frozen   |     |    | ▪              |                    |
| Immunohistology - Paraffin |     |    | ▪              |                    |
| ELISA                      | ▪   |    |                | 1/100 - 1/1000     |
| Western Blotting           |     |    | ▪              |                    |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

|                                       |   |
|---------------------------------------|---|
| <b>Target Species</b>                 | Human   |
| <b>Product Form</b>                   | Purified IgG - liquid   |
| <b>Preparation</b>                    | Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant |
| <b>Buffer Solution</b>                | Phosphate buffered saline   |
| <b>Preservative Stabilisers</b>       | <0.1% Sodium Azide (NaN <sub>3</sub> )  |
| <b>Approx. Protein Concentrations</b> | IgG concentration 1.0 mg/ml   |
| <b>Immunogen</b>                      | Purified light chains from human Ig.  |

|                                      |   |
|--------------------------------------|---|
| <b>External Database Links</b>       | <b>UniProt:</b><br><a href="#">P01834</a> <a href="#">Related reagents</a>  |
|                                      | <b>Entrez Gene:</b><br><a href="#">3514</a> IGKC <a href="#">Related reagents</a>   |
| <b>RRID</b>                          | AB_915299   |
| <b>Specificity</b>                   | <b>Mouse anti Human kappa light chain antibody, clone 4G7</b> recognizes a conformation-independent, constantly expressed epitope of the human immunoglobulin kappa chain.  |
| <b>Storage</b>                       | <p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p> |
| <b>Guarantee</b>                     | 12 months from date of despatch   |
| <b>Health And Safety Information</b> | Material Safety Datasheet documentation #10040 available at:<br><a href="https://www.bio-rad-antibodies.com/SDS/5268-6010">https://www.bio-rad-antibodies.com/SDS/5268-6010</a><br>10040  |
| <b>Regulatory</b>                    | For research purposes only  |

## Related Products

### Recommended Secondary Antibodies

|   |   |
|---|---|
| Goat Anti Mouse IgG (STAR77...)         | <a href="#">HRP</a>   |
| Rabbit Anti Mouse IgG (STAR12...)       | <a href="#">RPE</a>   |
| Goat Anti Mouse IgG (STAR70...)         | <a href="#">FITC</a>  |
| Goat Anti Mouse IgG IgA IgM (STAR87...) | <a href="#">Alk. Phos.</a> , <a href="#">HRP</a>  |
| Rabbit Anti Mouse IgG (STAR9...)        | <a href="#">FITC</a>  |
| Goat Anti Mouse IgG (STAR76...)         | <a href="#">RPE</a>   |
| Goat Anti Mouse IgG (H/L) (STAR117...)  | <a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> ,<br><a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> ,<br><a href="#">FITC</a> , <a href="#">HRP</a> |
| Rabbit Anti Mouse IgG (STAR13...)       | <a href="#">HRP</a>   |
| Goat Anti Mouse IgG (Fc) (STAR120...)   | <a href="#">FITC</a> , <a href="#">HRP</a>  |

### Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

|                                  |   |                  |   |               |   |
|----------------------------------|---|------------------|---|---------------|---|
| <b>North &amp; South America</b> | Tel: +1 800 265 7376<br>Fax: +1 919 878 3751<br>Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a> | <b>Worldwide</b> | Tel: +44 (0)1865 852 700<br>Fax: +44 (0)1865 852 739<br>Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a> | <b>Europe</b> | Tel: +49 (0) 89 8090 95 21<br>Fax: +49 (0) 89 8090 95 50<br>Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a> |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M410538:221028'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)