

Datasheet: 5276-2104

Description:	GOAT ANTI MOUSE IgM
Specificity:	IgM
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
ELISA	▪			1/20,000 - 1/100,000
Western Blotting	▪			1/2000 - 1/10,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

This antibody is suitable for use in a number of applications including immunodiffusion and immunoprecipitation.

Target Species	Mouse
Product Form	Purified Ig - liquid
Preparation	Purified Ig prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgM concentration 2.0 mg/ml
Immunogen	Mouse IgM mu heavy chain

**External Database
Links**

UniProt:

[P01873](#) [Related reagents](#)
[P01872](#) [Related reagents](#)

Entrez Gene:

[16019](#) Ighm [Related reagents](#)
[16019](#) Ighm [Related reagents](#)

RRID AB_620428

Specificity **Goat anti Mouse IgM polyclonal antibody** is specific for mouse IgM. Immunoelectrophoresis resulted in a single precipitin arc against anti-goat serum, mouse IgM and mouse serum. No reaction was observed against other mouse or human heavy or light chain proteins.

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 18 months from date of despatch.

**Health And Safety
Information** Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**North & South
America** Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M352444:190401'

Printed on 09 Feb 2021