

Datasheet: 5276-6504

BATCH NUMBER 162300

Description:	NATIVE PIG IgM
Name:	IgM
Format:	Purified
Product Type:	Antigen
Isotype:	IgM
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified porcine immunoglobulin from serum - liquid
Preparation	Delipidation, selective precipitation, molecular sieve chromatography
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	1.0 mg/ml
RRID	AB_615620

Product Information Single arc against anti-porcine IgM and porcine serum in IEP. No reaction was observed against anti porcine IgG-Fc. Some light chain cross reactivity will occur with anti-porcine IgG.

Purity	SDS PAGE: >95%
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10332 available at: https://www.bio-rad-antibodies.com/SDS/5276-6504 10332
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363259:200528'

Printed on 25 Mar 2023