

Datasheet: 5276-6504

BATCH NUMBER 165229

| | |
|----------------------|----------------|
| Description: | NATIVE PIG IgM |
| Name: | IgM |
| Format: | Purified |
| Product Type: | Antigen |
| Isotype: | IgM |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|---|
| Target Species | Pig |
| Product Form | Purified porcine immunoglobulin from serum - liquid |
| Preparation | Delipidation, selective precipitation, molecular sieve chromatography |
| Buffer Solution | TRIS buffered saline |
| Preservative Stabilisers | 0.1% Sodium Azide (NaN ₃) |
| Approx. Protein Concentrations | 1.0 mg/ml |
| RRID | AB_615620 |

Product Information Single arc against anti-porcine IgM and porcine serum in IEP. No reaction was observed against anti porcine IgG-Fc. Some light chain cross reactivity will occur with anti-porcine IgG.

| | |
|--------------------------------------|--|
| Purity | SDS PAGE: >95% |
| Storage | Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use. |
| Guarantee | 12 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10332 available at: https://www.bio-rad-antibodies.com/SDS/5276-6504 10332 |
| Regulatory | For research purposes only |

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M363259:200528'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)