

Datasheet: 5278-5204

BATCH NUMBER 171542

Description:	RAT ANTI HUMAN IgM MU CHAIN
Specificity:	IgM MU CHAIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	HM-7
Isotype:	IgM
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgM - liquid
Preparation	Purified IgM prepared by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgM concentration 1.0 mg/ml
Immunogen	Purified human IgM.

**External Database
Links**

UniProt:

[P01871](#) [Related reagents](#)

Entrez Gene:

[3507](#) IGHM [Related reagents](#)

RRID AB_618526

Specificity **Rat anti Human IgM mu chain antibody, clone HM-7** reacts with the mu chain of human IgM.

This antibody has not been cross-adsorbed and may react with homologous peptides from related species.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

**Health And Safety
Information** Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/5278-5204>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rat IgM (302001...) [Biotin](#), [RPE](#)

Goat Anti Rat IgM (STAR116...) [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M382528:210513'

Printed on 29 Jan 2026