

Datasheet: 5731-3000

BATCH NUMBER 161443

Description:	MOUSE ANTI HUMAN LUTEINIZING HORMONE RELEASING HORMONE RECEPTOR
Specificity:	LUTEINIZING HORMONE RELEASING HORMONE RECEPTOR
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	A9E4
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen	▪			
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein	IgG concentration 1.0 mg/ml

Concentrations

Immunogen Synthetic peptide corresponding to the amino acids 1 - 29 of the human luteinizing hormone releasing hormone receptor extracellular domain, conjugated to Bovine Serum Albumin.

External Database

Links

UniProt:[P30968](#)[Related reagents](#)**Entrez Gene:**[2798](#)

GNRHR

[Related reagents](#)

Synonyms

GRHR

RRID

AB_620154

Specificity

Mouse anti Human Luteinizing Hormone Releasing Hormone Receptor antibody, clone A9E4 recognizes human luteinizing hormone releasing hormone receptor (human GnRH receptor) found in the anterior pituitary.

References

1. Karande, A.A. *et al.* (1995) Establishment of immunological probes to study human gonadotropin-releasing hormone receptors. [Mol Cell Endocrinol. 114 \(1-2\): 51-6.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/5731-3000>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)

[HRP](#)

Rabbit Anti Mouse IgG (STAR12...)

[RPE](#)

Goat Anti Mouse IgG (STAR70...)

[FITC](#)

Goat Anti Mouse IgG IgA IgM (STAR87...)

[Alk. Phos.](#), [HRP](#)

Goat Anti Mouse IgG (STAR76...)

[RPE](#)

Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M381924:210512'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)