

Datasheet: 610001

**BATCH NUMBER 166139**

|                      |                             |
|----------------------|-----------------------------|
| <b>Description:</b>  | GOAT ANTI CHICKEN IgG (H/L) |
| <b>Specificity:</b>  | IgG (H/L)                   |
| <b>Other names:</b>  | IgY                         |
| <b>Format:</b>       | Purified                    |
| <b>Product Type:</b> | Polyclonal Antibody         |
| <b>Isotype:</b>      | Polyclonal IgG              |
| <b>Quantity:</b>     | 1 mg                        |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                     | Yes | No | Not Determined | Suggested Dilution |
|---------------------|-----|----|----------------|--------------------|
| ELISA               | ▪   |    |                |                    |
| Immunoprecipitation |     |    | ▪              |                    |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

|                                       |  |
|---------------------------------------|--|
| <b>Target Species</b>                 | Chicken  |
| <b>Product Form</b>                   | Purified IgG - liquid                            |
| <b>Preparation</b>                    | Purified IgG prepared by affinity chromatography |
| <b>Buffer Solution</b>                | Borate buffered saline                           |
| <b>Preservative Stabilisers</b>       | None present                                     |
| <b>Approx. Protein Concentrations</b> | IgG concentration 1.0mg/ml                       |
| <b>RRID</b>                           | AB_618828  |

**Specificity** **Goat anti Chicken IgG antibody** recognizes both the heavy and light chains of Chicken

IgG.

---

|                                      |  |
|--------------------------------------|--|
| <b>Storage</b>                       | Store at -20°C only.<br>Storage in frost-free freezers is not recommended.<br>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. |
| <b>Guarantee</b>                     | Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.   |
| <b>Health And Safety Information</b> | Material Safety Datasheet documentation #10123 available at:<br><a href="https://www.bio-rad-antibodies.com/SDS/610001">https://www.bio-rad-antibodies.com/SDS/610001</a><br>10123   |
| <b>Regulatory</b>                    | For research purposes only   |

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

|                                  |   |                  |   |               |   |
|----------------------------------|---|------------------|---|---------------|---|
| <b>North &amp; South America</b> | Tel: +1 800 265 7376<br>Fax: +1 919 878 3751<br>Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a> | <b>Worldwide</b> | Tel: +44 (0)1865 852 700<br>Fax: +44 (0)1865 852 739<br>Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a> | <b>Europe</b> | Tel: +49 (0) 89 8090 95 21<br>Fax: +49 (0) 89 8090 95 50<br>Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a> |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M401843:220718'

Printed on 01 Mar 2024