

Datasheet: 613005

BATCH NUMBER 163262

Description:	RABBIT F(ab') ₂ ANTI RAT IgG:HRP
Specificity:	IgG
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/50 - 1/200
ELISA	▪			1/4000 - 1/8000
Immunoblotting	▪			1/1000 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Rat

Product Form F(ab')₂ fragment of IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to rat IgG was raised by repeated immunisations of rabbits with rat IgG paraproteins. The preparation was purified by affinity chromatography on rat IgG covalently linked to agarose, and gel filtration chromatography of pepsin digested antibody.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 50% Glycerol

Immunogen Rat IgG.

