

Datasheet: 618501

BATCH NUMBER 166122

Description:	RABBIT ANTI RAT IgG (H/L)
Specificity:	IgG (H/L)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to rat IgG was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Rat IgG.

**External Database
Links**

UniProt:

P01835	Related reagents
P20766	Related reagents
P20760	Related reagents
P20759	Related reagents
P20761	Related reagents
P01836	Related reagents
P20767	Related reagents
P20762	Related reagents

Entrez Gene:

500180	LOC500180	Related reagents
363828	RGD1564318	Related reagents
299354	Ighg	Related reagents
679045	LOC679045	Related reagents
362795	LOC362795	Related reagents

RRID AB_618774

Specificity **Rabbit anti Rat IgG (H&L) antibody** recognizes the heavy and light chains of rat IgG1, IgG2a, IgG2b, IgG2c, and the light chains of rat IgM as demonstrated by ELISA.

Rabbit anti Rat IgG (H&L) antibody has been cross-absorbed with human sera and purified human paraproteins. It has minimal cross reactivity with human Immunoglobulins.

References 1. Potter, K.A. *et al.* (2012) Stab injury and device implantation within the brain results in inversely multiphasic neuroinflammatory and neurodegenerative responses. [J Neural Eng. 9: 046020.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10123 available at: <https://www.bio-rad-antibodies.com/SDS/618501>
10123

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M401832:220718'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)