

Datasheet: 6220-1004

BATCH NUMBER 163694

Description:	NATIVE HUMAN ALPHA 1 MICROGLOBULIN
Name:	ALPHA 1 MICROGLOBULIN
Other names:	AMBP
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Reconstitution	Reconstitute with 1.0ml sterile phosphate buffered saline (PBS). Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified Alpha 1 Microglobulin prepared by a combination of ion exchange, affinity and size exclusion chromatography.
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 1.0mg/ml after reconstitution.

**External Database
Links**

UniProt:

[P02760](#) [Related reagents](#)

Entrez Gene:

[259](#) AMBP [Related reagents](#)

Synonyms

HCP, ITIL

Product Information

Native Human alpha 1 microglobulin is native human Alpha-1-Microglobulin, a secreted peptide thought to play a role in the regulation of inflammatory processes.

This product runs at approximately 27 kDa under reducing conditions and as a dimer at approximately 54 kDa under non-reducing conditions.

Purity

>96% by SDS PAGE

References

1. Weyer, K. *et al.* (2013) Renal uptake of 99mTc-dimercaptosuccinic acid is dependent on normal proximal tubule receptor-mediated endocytosis. [J Nucl Med. 54 \(1\): 159-65.](#)

Storage

Prior to reconstitution store at +4°C.

After reconstitution store at -20°C.

Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed until date of expiry. Please see product label.

**Health And Safety
Information**

Material Safety Datasheet documentation #10268 available at:

<https://www.bio-rad-antibodies.com/SDS/6220-1004>

10268

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory

For research purposes only

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M412112:221109'

Printed on 19 Jan 2024