

Datasheet: 6220-1004

**BATCH NUMBER 163694**

<b>Description:</b>	NATIVE HUMAN ALPHA 1 MICROGLOBULIN
<b>Name:</b>	ALPHA 1 MICROGLOBULIN
<b>Other names:</b>	AMBP
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
<b>Reconstitution</b>	Reconstitute with 1.0ml sterile phosphate buffered saline (PBS). Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Purified Alpha 1 Microglobulin prepared by a combination of ion exchange, affinity and size exclusion chromatography.
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Total protein concentration 1.0mg/ml after reconstitution.

**External Database  
Links**

**UniProt:**

[P02760](#)    [Related reagents](#)

**Entrez Gene:**

[259](#)    AMBP    [Related reagents](#)

---

**Synonyms**

HCP, ITIL

---

**Product Information**

**Native Human alpha 1 microglobulin** is native human Alpha-1-Microglobulin, a secreted peptide thought to play a role in the regulation of inflammatory processes.

This product runs at approximately 27 kDa under reducing conditions and as a dimer at approximately 54 kDa under non-reducing conditions.

---

**Purity**

>96% by SDS PAGE

---

**References**

1. Weyer, K. *et al.* (2013) Renal uptake of 99mTc-dimercaptosuccinic acid is dependent on normal proximal tubule receptor-mediated endocytosis. [J Nucl Med. 54 \(1\): 159-65.](#)

---

**Storage**

Prior to reconstitution store at +4°C.

After reconstitution store at -20°C.

Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

Guaranteed until date of expiry. Please see product label.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10268 available at:  
<https://www.bio-rad-antibodies.com/SDS/6220-1004>  
10268

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

---

**Regulatory**

For research purposes only

---

**North & South  
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M412112:221109'

**Printed on 12 Aug 2023**