

Datasheet: 6220-1004

**BATCH NUMBER 172370**

<b>Description:</b>	NATIVE HUMAN ALPHA 1 MICROGLOBULIN
<b>Name:</b>	ALPHA 1 MICROGLOBULIN
<b>Other names:</b>	AMBP
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
<b>Reconstitution</b>	Reconstitute with 1.0ml sterile phosphate buffered saline (PBS). Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Total protein concentration 1.0mg/ml after reconstitution.
<b>External Database Links</b>	<b>UniProt:</b> <a href="https://www.uniprot.org/uniprot/P02760">P02760</a> <a href="#">Related reagents</a>

**Entrez Gene:**

[259](#) AMBP [Related reagents](#)

---

<b>Synonyms</b>	HCP, ITIL
-----------------	-----------

---

<b>Product Information</b>	<p><b>Native Human alpha 1 microglobulin</b> is native human Alpha-1-Microglobulin, a secreted peptide thought to play a role in the regulation of inflammatory processes.</p> <p>This product runs at approximately 27 kDa under reducing conditions and as a dimer at approximately 54 kDa under non-reducing conditions.</p>
----------------------------	---

---

<b>Purity</b>	>96% by SDS PAGE
---------------	------------------

---

<b>References</b>	1. Weyer, K. <i>et al.</i> (2013) Renal uptake of 99mTc-dimercaptosuccinic acid is dependent on normal proximal tubule receptor-mediated endocytosis. <a href="#">J Nucl Med. 54 (1): 159-65.</a>
-------------------	---

---

<b>Storage</b>	<p>This product is shipped at ambient temperature.</p> <p>Prior to reconstitution store at +4°C.</p> <p>After reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
------------------	--

---

<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10268 available at: <a href="https://www.bio-rad-antibodies.com/SDS/6220-1004">https://www.bio-rad-antibodies.com/SDS/6220-1004</a></p> <p>Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.</p> <p>As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious</p>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M438881:250523'

Printed on 29 Jan 2026