

Datasheet: 6240-0354

BATCH NUMBER 171414

Description:	MOUSE ANTI HUMAN BETA 2 MICROGLOBULIN
Specificity:	BETA 2 MICROGLOBULIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	B2M-01
Isotype:	IgG2a
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting		▪		

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Does not react with: Mouse, Dog, Rabbit, Chicken, Bovine
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	<0.1% Sodium Azide (NaN ₃)

Stabilisers

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Immunogen Beta-2 microglobulin from human urine.

External Database Links

UniProt:

[P61769](#) [Related reagents](#)

Entrez Gene:

[567](#) B2M [Related reagents](#)

RRID AB_617198

Specificity **Mouse anti Human beta 2 Microglobulin antibody, clone B2M-01** recognizes native beta 2 microglobulin which is a component of MHC class I molecules, which are present on all nucleated cells. Levels of beta 2 microglobulin can be used to assess renal function in kidney transplant recipients. Also known as B2M, it has been suggested that this protein can be used as a prognostic marker for HIV infection.

References 1. Hilgert, I. *et al.* (1984) The Use of Murine Monoclonal Antibody B2M-01 for Detection and Purification of Human beta 2- microglobulin. [Folia Biologica \(Praha\) 30: 369-376.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/6240-0354>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M395651:220519'

Printed on 29 Jan 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)