

Datasheet: 6400-2504

Description:	NATIVE MYCOPLASMA PNEUMONIAE
Name:	MYCOPLASMA PNEUMONIAE
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Inactivated <i>Mycoplasma pneumoniae</i> - liquid
Preparation	<i>Mycoplasma pneumoniae</i> cultured in a broth-based medium. Purified by washing, differential centrifugation and ether extraction. Residual solvent was removed by dry nitrogen. 6400-2504 has been inactivated by ether, but should be handled as if infectious.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Current, batch-specific concentration 2.1 mg/ml
Product Information	Native <i>Mycoplasma pneumoniae</i> is an antigen preparation containing soluble mycoplasma components in phosphate buffered saline.
Activity	Antigenic activity is 180% of internal reference standard.

Instructions For Use Sonicate immediately before use to ensure the preparation is uniform. This product may be used in a variety of immunoassay formats or may be further purified to meet the requirements of a particular assay format.

Storage Store at -70°C.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10209 available at:
<https://www.bio-rad-antibodies.com/SDS/6400-2504>
10209

The preparation has been inactivated by ether but should still be treated as infectious and handled appropriately

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M438215:250414'

Printed on 14 Apr 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)