

Datasheet: 640001

BATCH NUMBER 156324

Description:	SHEEP ANTI FITC
Specificity:	FITC
Other names:	FLUORESCEIN ISOTHIOCYANATE
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	
Immunofluorescence			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chemical
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography.
Antiserum Preparation	Antiserum to FITC was raised by repeated immunisation of sheep with highly purified antigen.
Buffer Solution	Borate buffered saline
Preservative	None present

Stabilisers

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Immunogen Fluorescein isothiocyanate (isomer 1).

RRID AB_620733

Specificity **Sheep anti FITC antibody** recognizes the fluorochrome Fluorescein Isothiocyanate (FITC). This reagent may be useful for amplification of staining using FITC conjugated reagents.

References

1. Heine, S. *et al.* (2011) CNGA3: A Target of Spinal Nitric Oxide/cGMP Signaling and Modulator of Inflammatory Pain Hypersensitivity. [J Neurosci. 31: 11184-92.](#)
2. Martínez-Sernández V *et al.* (2016) Usefulness of ELISA Methods for Assessing LPS Interactions with Proteins and Peptides. [PLoS One. 11 \(6\): e0156530.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10123 available at: <https://www.bio-rad-antibodies.com/SDS/640001>
10123

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350213:190307'

Printed on 19 Jan 2024