

Datasheet: 640005

**BATCH NUMBER 159261**

<b>Description:</b>	SHEEP ANTI FITC:HRP
<b>Specificity:</b>	FITC
<b>Other names:</b>	FLUORESCCEIN ISOTHIOCYANATE
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Chemical
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG was prepared from whole serum by affinity chromatography.
<b>Antiserum Preparation</b>	Antiserum to was raised by repeated immunisation of sheep with highly purified antigen.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Fluorescein isothiocyanate (isomer 1).

RRID AB\_619937

**Specificity** **Sheep anti FITC antibody** recognizes the fluorochrome Fluorescein Isothiocyanate (FITC). This reagent may be useful for amplification of staining using FITC conjugated reagents.

**References**

1. Heine, S. *et al.* (2011) CNGA3: A Target of Spinal Nitric Oxide/cGMP Signaling and Modulator of Inflammatory Pain Hypersensitivity. [J Neurosci. 31: 11184-92.](#)
2. Martínez-Sernández V *et al.* (2016) Usefulness of ELISA Methods for Assessing LPS Interactions with Proteins and Peptides. [PLoS One. 11 \(6\): e0156530.](#)
3. Aillaud, I. *et al.* (2022) A novel D-amino acid peptide with therapeutic potential (ISAD1) inhibits aggregation of neurotoxic disease-relevant mutant Tau and prevents Tau toxicity *in vitro.* [Alzheimers Res Ther. 14 \(1\): 15.](#)

**Storage** Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

**Guarantee** Guaranteed until date of expiry. Please see product label.

**Health And Safety Information** Material Safety Datasheet documentation #10106 available at: <https://www.bio-rad-antibodies.com/SDS/640005>  
10106

**Regulatory** For research purposes only

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350214:190307'

Printed on 01 Mar 2024