

Datasheet: 640005

BATCH NUMBER 163119

Description:	SHEEP ANTI FITC:HRP
Specificity:	FITC
Other names:	FLUORESCCEIN ISOTHIOCYANATE
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chemical
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG was prepared from whole serum by affinity chromatography.
Antiserum Preparation	Antiserum to was raised by repeated immunisation of sheep with highly purified antigen.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% Glycerol
Immunogen	Fluorescein isothiocyanate (isomer 1).

Specificity	Sheep anti FITC antibody recognizes the fluorochrome Fluorescein Isothiocyanate (FITC). This reagent may be useful for amplification of staining using FITC conjugated reagents.
References	<ol style="list-style-type: none"> 1. Heine, S. <i>et al.</i> (2011) CNGA3: A Target of Spinal Nitric Oxide/cGMP Signaling and Modulator of Inflammatory Pain Hypersensitivity. J Neurosci. 31: 11184-92. 2. Martínez-Sernández V <i>et al.</i> (2016) Usefulness of ELISA Methods for Assessing LPS Interactions with Proteins and Peptides. PLoS One. 11 (6): e0156530. 3. Aillaud, I. <i>et al.</i> (2022) A novel D-amino acid peptide with therapeutic potential (ISAD1) inhibits aggregation of neurotoxic disease-relevant mutant Tau and prevents Tau toxicity <i>in vitro.</i> Alzheimers Res Ther. 14 (1): 15. 4. Kenzel, J. <i>et al.</i> (2022) Selection of <i>Listeria monocytogenes</i> InIA-Binding Peptides Using Phage Display—Novel Compounds for Diagnostic Applications? Applied Microbiology. 2 (4): 921-933.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10106 available at: https://www.bio-rad-antibodies.com/SDS/640005 10106
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402092:220718'

Printed on 01 Mar 2024

