

Datasheet: 642005

BATCH NUMBER 155487

Description:	DONKEY ANTI GOAT IgG (H/L):HRP (RAT/MOUSE ADSORBED)
Specificity:	IgG (H/L)
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Goat
Species Cross Reactivity	Does not react with: Mouse, Rat
Product Form	Ig Fraction conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Donkey anti Goat IgG was prepared from pooled antisera from donkeys hyperimmunised with goat IgG. The antiserum was cross adsorbed against mouse and rat serum proteins to minimize cross-reactivity. The preparation was purified by affinity chromatography on goat IgG covalently linked to agarose.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% glycerol.

Immunogen	Goat IgG.
RRID	AB_609688
Specificity	Donkey anti Goat IgG antibody recognizes the heavy and light chains of goat IgG as demonstrated by ELISA. The preparation has minimal cross reaction with rat and mouse serum proteins, however, it may react with the light chains of other goat immunoglobulins.
References	<ol style="list-style-type: none"> Gerry, A.B. and Leake, D.S. (2014) Effect of low extracellular pH on NF-κB activation in macrophages. Atherosclerosis. 233: 537-44. Armour, K.L. <i>et al.</i> (2014) Clearance of Human IgG1-Sensitised Red Blood Cells In Vivo in Humans Relates to the In Vitro Properties of Antibodies from Alternative Cell Lines. PLoS One. 9: e109463.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10106 available at: https://www.bio-rad-antibodies.com/SDS/642005 10106
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350218:190307'

Printed on 01 Mar 2024