

Datasheet: 643002

BATCH NUMBER 169063

Description:	DONKEY ANTI RAT IgG (H/L):FITC (MOUSE ADSORBED)
Specificity:	IgG (H/L)
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat		
Species Cross Reactivity	Does not react with:Mouse		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)		
Approx. Protein Concentrations	1.0 mg/ml		
Immunogen	Rat IgG.		

**External Database
Links**

UniProt:

P20767	Related reagents
P20760	Related reagents
P20759	Related reagents
P20761	Related reagents
P20762	Related reagents
P01835	Related reagents
P01836	Related reagents
P20766	Related reagents

Entrez Gene:

679045	LOC679045	Related reagents
299354	Ighg	Related reagents
362795	LOC362795	Related reagents
500180	LOC500180	Related reagents
363828	RGD1564318	Related reagents

RRID AB_616788

Specificity **Donkey anti Rat IgG antibody** recognizes rat IgG heavy and light chains and demonstrates minimal cross reactivity with mouse serum proteins.

References

1. Bombardieri, M. *et al.* (2011) A BAFF/APRIL-dependent TLR3-stimulated pathway enhances the capacity of rheumatoid synovial fibroblasts to induce AID expression and Ig class-switching in B cells. [Ann Rheum Dis. 70 \(10\): 1857-65.](#)
2. Ivanescu, A.A. *et al.* (2015) Modifying Choroidal Neovascularization Development with a Nutritional Supplement in Mice. [Nutrients. 7 \(7\): 5423-42.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/643002>
10040

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

To

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

find a

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M428007:240301'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)