

Datasheet: 644005

BATCH NUMBER 156734

| | |
|----------------------|---|
| Description: | DONKEY ANTI RABBIT IgG (H/L):HRP (HUMAN/RAT/MOUSE ADSORBED) |
| Specificity: | IgG (H/L) |
| Format: | HRP |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | ▪ | | | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/4000 - 1/8000 |
| Western Blotting | | | ▪ | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Rabbit |
| Product Form | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid |
| Preparation | Purified IgG prepared by affinity chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 50% Glycerol |
| Immunogen | Rabbit IgG. |
| External Database Links | UniProt: |

[P01840](#) [Related reagents](#)
[P01839](#) [Related reagents](#)
[P01838](#) [Related reagents](#)
[P01847](#) [Related reagents](#)
[P01870](#) [Related reagents](#)
[P01841](#) [Related reagents](#)

Entrez Gene:

[100422811](#) IGKC2 [Related reagents](#)
[100009097](#) LOC100009097 [Related reagents](#)

RRID AB_619906

Specificity **Donkey anti rabbit IgG (H+L) antibody** recognizes rabbit immunoglobulin G heavy and light chains, and shows minimal cross-reactivity with human, rat and mouse serum proteins.
Donkey anti rabbit IgG (H+L) antibody may also recognize the light chains of other rabbit immunoglobulins.

References

1. Vega-Flores, G. *et al.* (2014) The GABAergic septohippocampal pathway is directly involved in internal processes related to operant reward learning. [Cereb Cortex. 24 \(8\): 2093-107.](#)
2. Belostotskaya, G. *et al.* (2015) Identification of cardiac stem cells within mature cardiac myocytes. [Cell Cycle. 17 : 1-8.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10106 available at:
<https://www.bio-rad-antibodies.com/SDS/644005>
10106

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)
[TMB CORE+ \(BUF062A\)](#)
[TMB SIGNAL+ \(BUF054A\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M350226:190307'

Printed on 01 Mar 2024