

Datasheet: 644005

**BATCH NUMBER 172784**

<b>Description:</b>	DONKEY ANTI RABBIT IgG (H/L):HRP (HUMAN/RAT/MOUSE ADSORBED)
<b>Specificity:</b>	IgG (H/L)
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Rabbit IgG.
<b>External Database Links</b>	<b>UniProt:</b>

[P01840](#) [Related reagents](#)  
[P01839](#) [Related reagents](#)  
[P01838](#) [Related reagents](#)  
[P01847](#) [Related reagents](#)  
[P01870](#) [Related reagents](#)  
[P01841](#) [Related reagents](#)

**Entrez Gene:**

[100422811](#) IGKC2 [Related reagents](#)  
[100009097](#) LOC100009097 [Related reagents](#)

---

**RRID** AB\_619906

---

**Specificity** **Donkey anti rabbit IgG (H+L) antibody** recognizes rabbit immunoglobulin G heavy and light chains, and shows minimal cross-reactivity with human, rat and mouse serum proteins.  
Donkey anti rabbit IgG (H+L) antibody may also recognize the light chains of other rabbit immunoglobulins.

---

**References**

1. Vega-Flores, G. *et al.* (2014) The GABAergic septohippocampal pathway is directly involved in internal processes related to operant reward learning. [Cereb Cortex. 24 \(8\): 2093-107.](#)
2. Belostotskaya, G. *et al.* (2015) Identification of cardiac stem cells within mature cardiac myocytes. [Cell Cycle. 17 : 1-8.](#)

---

**Storage** Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

---

**Health And Safety Information** Material Safety Datasheet documentation #10106 available at: <https://www.bio-rad-antibodies.com/SDS/644005>

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Useful Reagents

[TMB CORE \(BUF056A\)](#)  
[TMB CORE+ \(BUF062A\)](#)  
[TMB SIGNAL+ \(BUF054A\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M428009:240301'

**Printed on 29 Jan 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)