

Datasheet: 7750-0502

BATCH NUMBER 155613

Description:	SYNTHETIC HUMAN C-PEPTIDE
Name:	C-PEPTIDE
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified synthetic peptide - lyophilized
Reconstitution	Reconstitute with 0.5 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Solid phase synthesis
Buffer Solution	Acetate salt
Preservative Stabilisers	None present

External Database Links

UniProt:

[P01308](#)

[Related reagents](#)

Entrez Gene:

Product Information **Synthetic C-peptide** is a synthetic peptide corresponding to amino acids 3-33 of Human C-peptide, also known as 'Connecting peptide'.

C-peptide is the product of protease cleavage of proinsulin to produce insulin. Both Insulin and C-peptide are biologically active.

Protein Molecular Weight 3.02 kDa

Purity HPLC: >95%

Amino Acid Sequence EAEDLQVGQVELGGGPGAGSLQPLALEGSLQ

Storage Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 6 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10165 available at:
<https://www.bio-rad-antibodies.com/SDS/7750-0502>
10165

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M386893:210601'

Printed on 25 Mar 2023