

Datasheet: 7820-0604

BATCH NUMBER 158185

Description:	NATIVE HUMAN PROSTATE SPECIFIC ANTIGEN
Name:	PROSTATE SPECIFIC ANTIGEN
Other names:	KALLIKREIN 3
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human seminal fluid - liquid
Preparation	Sterile filtered through a 0.2 µM membrane
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	1.0 mg/ml

External Database Links

UniProt:

[P07288](https://www.uniprot.org/uniprot/P07288)

[Related reagents](#)

Entrez Gene:

Synonyms	APS
-----------------	-----

Product Information	<p>Native Human prostate specific antigen preparation is a purified form of human prostate-specific antigen (PSA), a glycoprotein of ~30 kDa found mainly in prostatic tissue and seminal fluid. PSA is also present in small quantities in the serum of normal men and is often elevated in the presence of prostate cancer and in other prostatic disorders.</p> <p>Native Human prostate specific antigen demonstrates high specific activity and low levels of any contaminating proteins, including those with trypsin-like protease activity (Manning et al. 2012).</p>
----------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Protein Molecular Weight	30 kDa
---------------------------------	--------

Purity	>40% by SDS-PAGE
---------------	------------------

References	<ol style="list-style-type: none">1. Manning, M,L, <i>et al.</i> (2012) Trypsin-like proteolytic contamination of commercially available psa purified from human seminal fluid. Prostate. 72 (11): 1233-8.2. Arya, S.K. and Bhansali, S. (2012) Anti-Prostate Specific Antigen (Anti-PSA) Modified Interdigitated Microelectrode-Based Impedimetric Biosensor for PSA Detection Biosensors Journal 1: H1106013. Chornokur, G. <i>et al.</i> (2011) Impedance-Based Miniaturized Biosensor for Ultrasensitive and Fast Prostate-Specific Antigen Detection Journal of Sensors, vol. 2011, Article ID 9837524. Manning ML <i>et al.</i> (2013) Proteolysis of complement factors iC3b and C5 by the serine protease prostate-specific antigen in prostatic fluid and seminal plasma. J Immunol. 190 (6): 2567-74.5. Piletska, E.V. <i>et al.</i> (2014) Microplates with enhanced immobilization capabilities controlled by a magnetic field J Chin Adv Mat Soc. 2: 118-29.6. Isaacs, J.T. <i>et al.</i> (2018) SYNTHETIC ENHANCEMENT OF THE T-CELL ARMAMENTARIUM AS AN ANTI-CANCER THERAPY US Patent application 20180148480 A17. Li, H. <i>et al.</i> (2020) Comprehensive role of prostate-specific antigen identified with proteomic analysis in prostate cancer J Cell Mol Med. 27 Jul. [Epub ahead of print].
-------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	<p>Material Safety Datasheet documentation #10258 available at: https://www.bio-rad-antibodies.com/SDS/7820-0604</p> <p>10258</p>
--------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M363449:200528'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)