

Datasheet: 7861-1079

Description:	NATIVE HUMAN PROTEIN S
Name:	PROTEIN S
Format:	Purified
Product Type:	Antigen
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human plasma - liquid
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	1mM Benzamidine
Approx. Protein Concentrations	1.0 mg/ml
Purity	SDS-PAGE: >95% Single band at 69 kDa on SDS-PAGE gels, and a doublet upon incubation with 2-mercaptoethanol.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10301 available at:
10301: <https://www.bio-rad-antibodies.com/uploads/MSDS/10301.pdf>

Donor material tested and found negative for HIV1/2 antibodies, HCV, Rapid Plasma Reagin and HBsAg.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M363456:200528'

Printed on 29 Aug 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)