

Datasheet: 7863-2004

Description:	MOUSE ANTI HUMAN PROTEIN GENE PRODUCT 9.5
Specificity:	PROTEIN GENE PRODUCT 9.5
Other names:	PGP 9.5, UCHL1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	13C4
Isotype:	IgG2a
Quantity:	0.2 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin (1)	▪			1/200
ELISA	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

(1) Antigen is stable in formalin fixed paraffin embedded sections, however we recommend fixation in 95% ethanol/5% acetic acid for 2-3 hours prior to paraffin embedding. Can be used without acid/ethanol fixation if the sections are subjected to microwave treatment in citrate buffer by standard methods.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Rat, Guinea Pig, Pig, Beaver, Donkey</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid

Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	1 mg/ml
Immunogen	Native, from brain
External Database Links	<p>UniProt: P09936 Related reagents</p> <p>Entrez Gene: 7345 UCHL1 Related reagents</p>
RRID	AB_620255
Specificity	<p>Mouse anti Human protein gene product 9.5 antibody, clone 13C4 recognizes protein gene product 9.5 (PGP9.5), a ubiquitin hydrolase widely expressed in neuronal tissues representing, by weight 1-2% of total soluble brain proteins. PGP9.5, also known as ubiquitin C-terminal hydrolase 1 (UCHL-1), is involved in the regulation of the ubiquitin pathway.</p> <p>Mouse anti Human protein gene product 9.5 antibody, clone 13C4 stains neuronal cell bodies and axons in the CNS and periphery, small nerve fibres in peripheral tissues, neuroendocrine cells in the pituitary, thyroid, pancreas and tumours of the DNES. Mouse anti Human protein gene product 9.5 antibody, clone 13C4 stains neuroendocrine cells in human adult gut (unlike Mouse anti Human Protein Gene Product 9.5 antibody, clone 31A3).</p> <p>Clones 31A3 and 13C4 each recognize a different epitope towards the N-terminus of the protein.</p>
References	<ol style="list-style-type: none"> 1. Wilson, P.O. <i>et al.</i> (1988) The immunolocalization of protein gene product 9.5 using rabbit polyclonal and mouse monoclonal antibodies. Br J Exp Pathol. 69 (1): 91-104. 2. Cooke, H.J. <i>et al.</i> (1999) Activation of neuronal adenosine A1 receptors suppresses secretory reflexes in the guinea pig colon. Am J Physiol. 276 (2 Pt 1): G451-62. 3. Sienkiewicz, W. <i>et al.</i> (2000) Has active immunization against gonadotrophin-releasing hormone any effect on testis innervation in the pig? An immunohistochemical study. Anat Histol Embryol. 29 (4): 247-54. 4. Sasaki, H. <i>et al.</i> (2001) Expression of the protein gene product 9.5, PGP9.5, is

- correlated with T-status in non-small cell lung cancer. [Jpn J Clin Oncol. 31: 532-5.](#)
5. Komori, N. *et al.* (2003) Presence of beta-arrestin-1 immunoreactivity in the cutaneous nerve fibers of rat glabrous skin. [Brain Res. 988 \(1-2\): 121-9.](#)
 6. Kaleczyc, J. *et al.* (2007) The distribution and chemical coding of intramural neurons supplying the porcine stomach - the study on normal pigs and on animals suffering from swine dysentery. [Anat Histol Embryol. 36 \(3\): 186-93.](#)
 7. Pidsudko, Z. *et al.* (2008) Distribution and chemical coding of intramural neurons in the porcine ileum during proliferative enteropathy. [J Comp Pathol. 138 \(1\): 23-31.](#)
 8. Kotani, T. *et al.* (2010) Expression of PTPRO in the interneurons of adult mouse olfactory bulb. [J Comp Neurol. 518: 119-36.](#)
 9. Godlewski J & Pidsudko Z (2012) Characteristic of galaninergic components of the enteric nervous system in the cancer invasion of human large intestine. [Ann Anat. 194 \(4\): 368-72.](#)
 10. Dudek, A. *et al.* (2012) Immunohistochemical characterization of neurons in the vestibular ganglion (Scarpa's ganglion) of the pig. [Pol J Vet Sci.15: 499-507.](#)
 11. Bulc, M. *et al.* (2013) Immunohistochemical characterization of the porcine nodose ganglion. [Acta Histochem. 115 \(5\): 440-6.](#)
 12. Burliński, P.J. (2012) Inflammation- and axotomy-induced changes in cocaine- and amphetamine-regulated transcript peptide-like immunoreactive (CART-LI) nervous structures in the porcine descending colon. [Pol J Vet Sci. 15 \(3\): 517-24.](#)
 13. Buels, K.S. *et al.* (2012) Non-bronchodilating mechanisms of tiotropium prevent airway hyperreactivity in a guinea-pig model of allergic asthma. [Br J Pharmacol. 165: 1501-14.](#)
 14. Akazawa, N. *et al.* (2014) Neuroendocrine carcinoma of the esophagus: clinicopathologic study of 10 cases and verification of the diagnostic utility of mASH1, NeuroD1, and PGP9.5 [Esophagus. 11 \(4\): 245-57.](#)
 15. Zalecki, M. (2015) The Influence of Antral Ulcers on Intramural Gastric Nerve Projections Supplying the Pyloric Sphincter in the Pig (*Sus scrofa domestica*)-Neuronal Tracing Studies. [PLoS One. 10 \(5\): e0126958.](#)
 16. Szymanska, K. *et al.* (2018) The Influence of High and Low Doses of Bisphenol A (BPA) on the Enteric Nervous System of the Porcine Ileum. [Int J Mol Sci. 2018;19 \(3\): 917.](#)
 17. Zalecki, M. *et al.* (2019) Enteric nervous system in the European beaver (*Castor fiber*) pylorus - an immunohistochemical study. [Pol J Vet Sci. 22 \(1\): 101-7.](#)
 18. Li, T.T. *et al.* (2020) Effects of different culture systems on the culture of prepuberal buffalo (*Bubalus bubalis*) spermatogonial stem cell-like cells *in vitro*. [J Vet Sci. 21 \(1\): e13.](#)
 19. Choi, Y. *et al.* (2020) Stage-Dependent Expression of Protein Gene Product 9.5 in Donkey Testes. [Animals \(Basel\). 10 \(11\): 2169.](#)
 20. Zalecki, M. *et al.* (2020) Inferior vagal ganglion galaninergic response to gastric ulcers. [PLoS One. 15 \(11\): e0242746.](#)
 21. Zalecki, M. *et al.* (2021) Cocaine and Amphetamine Regulated Transcript (CART) Expression Changes in the Stomach Wall Affected by Experimentally Induced Gastric Ulcerations. [Int J Mol Sci. 22 \(14\):7437.](#)
 22. Makowska, K. *et al.* (2022) The Comparison of the Influence of Bisphenol A (BPA) and Its Analogue Bisphenol S (BPS) on the Enteric Nervous System of the Distal Colon in Mice [Nutrients. 15 \(1\): 200.](#)
 23. Makowska, K. *et al.* (2023) Bisphenol A Effects on Neurons' Neurochemical Character in the Urinary Bladder Intramural Ganglia of Domestic Pigs. [Int J Mol Sci.24 \(23\): 16792.](#)

24. Rytel, L. *et al.* (2021) Chemically-Induced Inflammation Changes the Number of Nitroergic Nervous Structures in the Muscular Layer of the Porcine Descending Colon. [Animals \(Basel\). 11 \(2\): 394.](#)
25. Shakeel, M. *et al.* (2022) Seasonal changes in the expression of molecular markers of stallion germ cells. [J Equine Vet Sci. 118: 104109.](#)
26. Franke-Radowiecka, A. *et al.* (2020) Nerve structures of the heart and their immunohistochemical characterization in 10-week-old porcine fetuses. [C R Biol. 343 \(1\): 53-62.](#)
27. Bulc, M. *et al.* (2022) Administration of Different Doses of Acrylamide Changed the Chemical Coding of Enteric Neurons in the Jejunum in Gilts. [Int J Environ Res Public Health. 19 \(21\):14514.](#)
28. Makowska, K. *et al.* (2023) Effects of the long-term influence of bisphenol A and bisphenol S on the population of nitroergic neurons in the enteric nervous system of the mouse stomach. [Sci Rep. 13 \(1\): 331.](#)
29. Kawabata, Y. *et al.* (2024) Hair Follicle Damage after 100 mGy Low-Dose Fractionated X-Ray Irradiation and the Protective Effects of TEMPOL, a Stable Nitroxide Radical, against Radiation. [Radiat Res. 201 \(2\): 115-25.](#)
30. Makowska, K. *et al.* (2021) Bisphenol A affects vipergic nervous structures in the porcine urinary bladder trigone. [Sci Rep. 11 \(1\): 12147.](#)
31. Rychlik, A. *et al.* (2020) The T2 Toxin Produced by *Fusarium* spp. Impacts Porcine Duodenal Nitric Oxide Synthase (nNOS)-Positive Nervous Structures-The Preliminary Study. [Int J Mol Sci. 21 \(14\): 5118.](#)
32. Makowska, K. *et al.* (2023) Influence of bisphenol A and its analog bisphenol S on cocaine- and amphetamine-regulated transcript peptide-positive enteric neurons in the mouse gastrointestinal tract. [Front Mol Neurosci. 16: 1234841.](#)
33. Brzozowska, M. & Całka, J. (2023) Acetylsalicylic Acid Supplementation Affects the Neurochemical Phenotyping of Porcine Duodenal Neurons. [Int J Mol Sci. 24 \(12\) :9871.](#)
34. Baert, Y. *et al.* (2020) A multi-organ-chip co-culture of liver and testis equivalents: a first step toward a systemic male reprotoxicity model. [Hum Reprod. 35 \(5\): 1029-44.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/7863-2004>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M405301:220916'

Printed on 12 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)