

Datasheet: 7895-0104

BATCH NUMBER 170411

Description:	NATIVE RABBIT PYRUVATE KINASE
Name:	PYRUVATE KINASE
Format:	Purified
Product Type:	Purified Protein
Quantity:	50 KU

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Semi purified enzyme from rabbit muscle - lyophilised
Reconstitution	Reconstitute with sterile distilled water. For long term storage, addition of a preservative is recommended.
Preservative Stabilisers	None present
Approx. Protein Concentrations	75%

External Database Links

UniProt:

[P11974](#)

[Related reagents](#)

Entrez Gene:

[100008676](#)

LOC100008676

[Related reagents](#)

Activity 375 U/mg solid
500 U/mg protein

1 unit liberates 1 micromole pyruvate per minute from PEP in the presence of ADP at 30°C and pH7.6.

Storage Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Avoid repeated freezing and thawing as this may denature the protein.
Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10268 available at:
<https://www.bio-rad-antibodies.com/SDS/7895-0104>
10268

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363465:200528'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)