

Datasheet: 7970-0504

**BATCH NUMBER 167575**

<b>Description:</b>	NATIVE HUMAN RETINOL BINDING PROTEIN
<b>Name:</b>	RETINOL BINDING PROTEIN
<b>Other names:</b>	RBP
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
<b>Reconstitution</b>	Reconstitute with 1.0 ml sterile phosphate buffered saline. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	1.0 mg/ml after reconstitution
<b>External Database Links</b>	<b>UniProt:</b> <a href="https://www.uniprot.org/entry/P02753">P02753</a> <a href="#">Related reagents</a>

**Entrez Gene:**

[5950](#) RBP4 [Related reagents](#)

---

**Product Information** **Native Human retinol binding protein** is a purified human retinol binding protein preparation. Retinol binding protein is a secreted peptide that delivers retinol from liver stores to the peripheral tissues.

---

**Protein Molecular Weight** 21 kDa

---

**Purity** >96% by SDS PAGE

---

**Storage** Prior to reconstitution store at +4°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein.

---

**Guarantee** Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information** Material Safety Datasheet documentation #10268 available at:  
<https://www.bio-rad-antibodies.com/SDS/7970-0504>  
10268

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M426272:231207'

**Printed on 10 Jul 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)