

Datasheet: 8129-9804

BATCH NUMBER 147167

Description:	GOAT ANTI ROTAVIRUS
Specificity:	ROTAVIRUS
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
Western Blotting			▪	
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Viral
Product Form	Ig fraction - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	4.0 mg/ml
Immunogen	Nebraska calf diarrhea virus.
RRID	AB_618846

Specificity

Goat anti Rotavirus antibody reacts with rotavirus and is specific for ICPs and late structural (virion) antigens. It has been reported to exhibit 90% cross reactivity with human rotavirus and does not react with HEp-2 cells or WI-38 cells in immunofluorescence.

Rotavirus is a virus that infects the bowels and is the most common cause of gastroenteritis in infants. 7 major groups have currently been identified, 3 of which infect humans (groups A, B and C). The rotavirus genome consists of 11 segments of double stranded RNA held in the inner core of the three layered, icosahedral virion particle. The genome codes for 6 virus proteins (VP1,2,3,4,6,7) and 6 non-structural proteins (NSP1-6). Outer capsid proteins are VP4 and VP7, inner capsid -VP6 and core -VP2. The immunity acquired from exposure to rotavirus appears to be type specific following initial infection; therefore, multiple serotypes of rotavirus mean multiple opportunities for infection.

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/8129-9804>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363476:200528'

Printed on 01 Mar 2024