

Datasheet: 8203-0050

BATCH NUMBER 220414

| | |
|----------------------|--------------------------------------|
| Description: | RABBIT ANTI SACCHAROMYCES CEREVISIAE |
| Specificity: | SACCHAROMYCES CEREVISIAE |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|--------------------------|-----|----|----------------|--------------------|
| Immunohistology - Frozen | | | ▪ | |
| ELISA | | | ▪ | |
| Western Blotting | ▪ | | | |
| Immunofluorescence | | | ▪ | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|--|
| Target Species | Yeast |
| Product Form | Purified IgG - liquid |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) |
| Approx. Protein Concentrations | 1 mg/ml |
| Immunogen | Whole cells of <i>S. cerevisiae</i> . |
| RRID | AB_620518 |

| | |
|--------------------------------------|--|
| Specificity | Rabbit anti <i>Saccharomyces cerevisiae</i> antibody recognises <i>Saccharomyces cerevisiae</i> , a rapidly growing budding yeast commonly used in brewing and baking. <i>S. cerevisiae</i> has been studied intensively as a eukaryotic cell model as it is easy to culture but maintains the internal complexity of plant and animal cells. |
| Storage | Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. |
| Guarantee | 12 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/8203-0050 10040 |
| Regulatory | For research purposes only |

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363479:200528'

Printed on 30 Jan 2024