

Datasheet: 8210-0407

BATCH NUMBER 157369

Description:	MOUSE ANTI SALMONELLA TYPHIMURIUM LPS
Specificity:	SALMONELLA TYPHIMURIUM LPS
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	1E6
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified LPS from <i>S. typhimurium</i> .
RRID	AB_619235
Specificity	<p>Mouse anti <i>Salmonella typhimurium</i> LPS antibody, clone 1E6 reacts with an unidentified LPS determinant of <i>Salmonella typhimurium</i>.</p> <p>Clone 1E6 has been shown to not cross react with <i>S. paratyphi A</i>, <i>S. choleraesuis</i>, <i>S. newport</i>, <i>S. enteritidis</i>, <i>S. anatum</i>, <i>S. selandia</i>, <i>E. coli</i> 055:B5, <i>E. coli</i> K12, <i>Klebsiella pneumoniae</i>.</p>
References	<p>1. Rakebrandt, N. <i>et al.</i> (2014) Antibody- and TRIM21-dependent intracellular restriction of <i>Salmonella enterica</i>. Pathog Dis. 72 (2): 131-7.</p> <p>2. Schenk, F. <i>et al.</i> (2018) Development of a paper-based lateral flow immunoassay for simultaneous detection of lipopolysaccharides of <i>Salmonella</i> serovars. Anal Bioanal Chem. 410 (3): 863-8.</p>
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/8210-0407</p> <p>10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M375001:201211'

Printed on 11 Jul 2024