

Datasheet: 8429-9004

BATCH NUMBER 170128

Description:	STREPTAVIDIN
Name:	STREPTAVIDIN
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Product Form	Purified streptavidin from <i>Streptomyces avidinii</i> - lyophilized
---------------------	-----------------------------------------------------------------------

Reconstitution	Reconstitute with distilled water or saline. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
-----------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparation	Purified streptavidin prepared by affinity chromatography
--------------------	-----------------------------------------------------------

Buffer Solution	0.01M Sodium Phosphate 0.01M Sodium Borate 0.15M Sodium Chloride
------------------------	------------------------------------------------------------------------

Preservative	1% Mannitol
Stabilisers	1% Dextran

Protein Molecular Weight	53 kDa
---------------------------------	--------

Activity	17 U/mg streptavidin One unit will bind 1 µg of D-biotin.
-----------------	--------------------------------------------------------------

Purity	>95% by SDS-PAGE
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #20485 available at: https://www.bio-rad-antibodies.com/SDS/8429-9004
Regulatory	For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M429173:240320'

Printed on 20 Feb 2026