

Datasheet: 8900-1004

BATCH NUMBER 159934

Description:	NATIVE HUMAN THYROGLOBULIN
Name:	THYROGLOBULIN
Other names:	Tg
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human thyroid glands - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.0 mg/ml
External Database Links	<p>UniProt: P01266 Related reagents</p> <p>Entrez Gene: 7038 TG Related reagents</p>

Protein Molecular Weight	660 kDa, confirmed by chromatography and electrophoresis (under certain conditions the 330 kDa subunit may be observed)
Purity	>96% by SDS PAGE
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10209 available at: https://www.bio-rad-antibodies.com/SDS/8900-1004 10209 Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies. As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M362551:200514'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)