

Datasheet: 8900-1354

BATCH NUMBER 156902

| | |
|----------------------|-----------------------------|
| Description: | NATIVE BOVINE THYROGLOBULIN |
| Name: | THYROGLOBULIN |
| Other names: | Tg |
| Format: | Purified |
| Product Type: | Purified Protein |
| Quantity: | 100 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------|-----|----|----------------|--------------------|
| ELISA | ■ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|---|
| Target Species | Bovine |
| Product Form | Purified protein from bovine thyroid glands - liquid |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | None present |
| Approx. Protein Concentrations | Total protein concentration 5.0 mg/ml |
| External Database Links | <p>UniProt: P01267 Related reagents</p> <p>Entrez Gene: 280706 TG Related reagents</p> |

Product Information Bovine thyroglobulin, a precursor of the iodated thyroid hormones thyroxine (T4) and triiodothyronine (T3) that affect protein, fat and carbohydrate metabolism. Thyroglobulin is produced by thyrocytes within the thyroid gland and is stored within the follicular lumen.

References 1. Uyttebroek, L. *et al.* (2013) Expression of neuropeptides and anoctamin 1 in the embryonic and adult zebrafish intestine, revealing neuronal subpopulations and ICC-like cells. [Cell Tissue Res. 354 \(2\): 355-70.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10209 available at: <https://www.bio-rad-antibodies.com/SDS/8900-1354>
10209

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363536:200528'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)